

TEN-MINUTE GENOMIC DNA PREPARATION FROM YEAST

Reagents:

- YPD
- Sterile water
- Lysis buffer: 2% Triton X-100, 1% SDS, 100 mM NaCl, 10mM Tris-HCl (pH 8.0),
1mM Na₂EDTA.
- Phenol:Chloroform:isoamyl alcohol (25:24:1).
- Glass beads: 0.45-0.5 mm beads (Sigma) soaked in nitric acid and washed in
distilled water. Beads should be dried before use.
- TE (pH 8.0): 10 mM Tris-HCl, 1 mM Na₂EDTA
- 100% EtOH
- Rnase A (10 mg/ml)
- 7.5 M ammonium acetate

Protocol:

1. Grow 10-ml yeast culture to saturation in YPD at 30°C.
2. Collect the cells by centrifugation for two minutes in a clinical centrifuge.
Remove the supernatant and resuspend the cells in 0.5 ml of distilled
water. Transfer the cells to a 1.5-ml microfuge tube and collect them by
centrifugation for 5 seconds in a microfuge,
3. Decant the supernatant and briefly vortex the tube to resuspend the pellet
in the residual liquid.
4. Add 0.2 ml of lysis buffer and 0.2 ml phenol:chloroform:isoamyl alcohol.
Then add 0.3 g of glass beads.

5. Vortex for 3-4 minutes. Then add 0.2 ml of TE (pH 8.0)
6. Centrifuge for 5 minutes in a microfuge. Transfer the aqueous layer to a fresh tube. Add 1 ml of 100% ethanol. Mix by inversion.
7. Centrifuge for 2 minutes in a microfuge. Discard the supernatant.
Resuspend the pellet in 0.4 ml of TE plus 3 μ l of a 10 mg/ml solution of RNase A. Incubate for 5 minutes at 37°C. Add 5.4 μ l of 7.5 M ammonium acetate plus 1 ml of 100% ethanol. Mix by inversion.
8. Centrifuge for 2 minutes in a microfuge. Discard the supernatant. Air-dry the pellet and resuspend in 50 μ l of TE.
9. Use 10 μ l for each sample to be analyzed by Southern blotting. This is ~2-4 μ g of DNA.