

## **FeBABE-mediated hydroxyl radical cleavage assay**

### **DNA hybridization**

1. DNA containing at least two phosphorothioates at alternate bases can be obtained from Invitrogen.

The following codes are used to modify bases

A-Phosphorothioate =F

C-Phosphorothioate =O

G-Phosphorothioate =E

T-Phosphorothioate =Z

2. All oligos mentioned here have the same bottom strand annealed to different top strands. The top strands contain biotin on the 5' end.
3. Anneal oligos in 100  $\mu$ l of 1X TE+ 0.1M NaCl by heating at 100°C for 5 mins and then slow cooling overnight. Add 10 $\mu$ g of biotinylated top strand and 20 $\mu$ g of the bottom strand so that it is always forced into the hybridized state
4. Ethanol precipitate the annealed oligo by adding 6  $\mu$ l of 3M NaOAc and 300  $\mu$ l of 100% EtOH to 100  $\mu$ l of annealed oligo. Store the pellet at -20°C for at least 4 hours and then wash with ~500  $\mu$ l of 80% EtOH and dry at room temperature.
5. Resuspend the pellet in 20mM MOPS pH 7.9 to a final concentration of 300 ng/  $\mu$ l.

### **Coupling to FeBABE**

6. FeBABE is obtained from Dojindo (catalog no: F279-10). Resuspend the yellowish powder in 50  $\mu$ l of water (to a final concentration of 34mM) and store at -80 °C till needed.
7. The coupling reaction is done in 20mM MOPS pH 7.9. To 600 ng of the oligo, 3.5  $\mu$ l of 34mM of FeBABE is added.

2  $\mu$ l hybridized biotinylated oligo

3.5  $\mu$ l of FeBABE (34mM)

0.8  $\mu$ l of 0.5M MOPS pH7.9

13.7  $\mu\text{l}$  of  $\text{H}_2\text{O}$

8. Incubate at  $50^\circ\text{C}$  for 16 hours (Put the reaction in the thermocycler and use the program INC under user ROS)

### **Coupling FeBABE conjugated DNA to beads**

9. Wash 40  $\mu\text{l}$  of Dynabeads M-280 Streptavidin (catalog no: 112.05D) for each reaction with 1X B&W buffer (manufacturer recommended) twice and once with 20mM MOPS pH 7.9
10. Pull down the beads on a magnet and remove the washes. Resuspend the beads in 300  $\mu\text{l}$  of 20mM MOPS pH7.9
11. Add the FeBABE-coupled DNA to the beads and let it incubate by rotating at  $4^\circ\text{C}$  for 4 hours
12. Then pull down the beads, remove the supernatant and wash the beads with 20mM MOPS pH7.9 twice to get rid of the unbound DNA and the excess FeBABE reagent.
13. Then, wash the beads once with the final reaction buffer (1X Reaction buffer -DTT ).

### **FEABE cleavage reaction**

14. Pull down and remove the wash. Resuspend the beads in 300  $\mu\text{l}$  of Reaction mix+water (containing DTT).
15. Add 6  $\mu\text{l}$  of 1  $\mu\text{M}$  TBP (20nM final concentration). Incubate the tubes on a rotator at room temperature for 20 minutes.
16. Pull down the beads using the magnet and remove the flow-through.

17. Resuspend the beads in 285  $\mu\text{l}$  of Reaction mix+water. Add 15  $\mu\text{l}$  of 145nM Mot1 and incubate at room temperature on the rotator for 30 minutes.
18. Then pull down the beads and remove the supernatant.
19. Wash the beads with 1 ml of 1X Reaction buffer (-DTT) once. (If reactions need to be split, this is a good point to split it.)
20. Pull down the beads on the magnet and remove the wash completely.
21. Resuspend the beads in 7.5  $\mu\text{L}$  of 1X Reaction buffer (-DTT). Add 2.5  $\mu\text{L}$  of 50% glycerol to each. Spin down quickly to bring down the beads (The glycerol is added to localize the free radicals; DTT is a radical scavenger, so you don't want it at this stage.)
22. To do the reaction, add 1.25  $\mu\text{L}$  of 50mM sodium ascorbate. Mix and then quickly add 1.25  $\mu\text{L}$  of 50mM  $\text{H}_2\text{O}_2$ /10mM EDTA. React for 5 minutes at 37°C (This can vary according to the cleavage required. I have tried 5 mins at 30°C or 37°C)
23. Quench with 6  $\mu\text{L}$  of 4X SB with 1  $\mu\text{L}$  of 1M DTT.
24. Separate the cleaved fragments using SDS-PAGE (I used 14% gels for Mot1). Western blot using an antibody that recognizes one terminus of the protein (For Mot1, I used the antibody that recognized the C terminus (Auble et al., 1997))

**1X binding and washing (B&W) buffer**

5mM Tris HCl pH7.5

0.5mM EDTA

1M NaCl

**5X Reaction buffer (same as 5X binding buffer)**

20% glycerol                      400 $\mu$ L of 50% glycerol

20mM Tris-Cl, pH 8              20 $\mu$ L of 1M

300mM KCl                        300 $\mu$ L of 1M

25mM MgCl<sub>2</sub>                      25  $\mu$ L of 1M

500mg/ml BSA                    10  $\mu$ L of 50mg/ml

245  $\mu$ L water

**1X Reaction buffer (for washes)**

Dilute 5X reaction buffer in water to 1X

**Reaction mix (For 300  $\mu$ L)**

60  $\mu$ L 5X binding buffer

3  $\mu$ L 100mM DTT

30  $\mu$ L 1% Brij

Adjust volume with water

**50mM H<sub>2</sub>O<sub>2</sub>/10mM EDTA**

974  $\mu$ L water

20  $\mu$ L 500mM EDTA

5.7  $\mu$ L of 30% H<sub>2</sub>O<sub>2</sub>

**50mM ascorbate**

99mg of sodium ascorbate in 10mL of water. (You can store the ascorbate at -80°C. But I just remake it every time as it is easier than to thaw it.)

**ADP-AIF<sub>4</sub>**

1  $\mu$ L of 2mM ADP

0.5  $\mu$ L of 50mM NaF

0.1  $\mu$ L of 10 mM AlCl<sub>3</sub>

Mix well and add 1.6  $\mu$ L to the reaction

**4X SB+DTT**

To 6  $\mu$ L of 4XSB, add 1  $\mu$ L of 1M DTT (Add that to quench the reaction)