

Northern Blots

(based on protocol from Mitch Smith Lab)

Use gel apparatus designated for RNA. Wipe apparatus with "RNase Away".

1% Agarose gel (100ml)

73 ml sterile water

1 g agarose

Boil in microwave and cool a little

Add 10 ml 10X MOPS and 17 ml 37% Formaldehyde

Pour.

10X MOPS

0.2 M MOPS 23.1 g/500ml

50 mM Sodium Acetate 8.3 ml of 3M/500 ml

10mM EDTA 10ml of 0.5M/500ml

adjust pH 7.0

Sample preparation and electrophoresis

Add 20 µg each RNA to eppendorf tube, dry samples in speed-vac. Add 17 µl Sample Buffer, heat samples at 65 C for 10 min, add 3 µl 10X dye, and then load gel and run in 1X MOPS

Run gel 1-2 hrs at 90 V

Photograph gel and then transfer

Sample Buffer

50% Formamide 500 µl/ml

2.2 M Formaldehyde (37%) 170 µl/ml

1X MOPS 100 µl of 10X/ml

10X Dye Solution

50% glycerol 100 µl

0.3% Bromophenol Blue 6 µl (10%)

Ethidium Bromide 25 µl 4mg/ml

Sterile water 75 µl

Transfer


Rinse gel twice for 15 min each in sterile water.

Cut nitrocellulose (Schleicher & Schuell Protran BA85, #10402580) to size of gel. Float nitrocellulose in water to wet and then soak in 20X SSC.

20X SSC

3.0M NaCl 175g/L
 0.3 M sodium citrate 88g/L
 Autoclave using liquid cycle for at least 40 minutes

Blot gel by capillary method using a reservoir of 20X SSC (see figure), usually overnight.



Place blot on plastic wrap, then UV crosslink RNA to blot using Stratalinker on automode.

Hybridization

Prehybridize blot at 42 C for 1-2 hrs in Hybridization Solution (plus freshly boiled herring sperm DNA (100 µg/ml final) and 1 ml of 50% dextran sulfate per 5 ml of Hybridization Solution) Use roller bottle in hybridization oven.

Hybridization Solution (100ml, store at 4 C)

50 ml formamide
 1 ml 10% SDS
 25 ml 20X SSC
 5 ml 1M sodium phosphate, pH 6.8
 10 ml 50X Denhardt's solution
 9 ml H2O

Before use:

Add boiled SS DNA to 0.1 mg/ml

Add 20 ml 50% dextran sulfate per 100 ml

50% Dextran sulfate: 25 g/50 ml water, mix overnight and store at 4 C

Fish Sperm DNA (10 mg/ml, store at -20C)

5030 μ l H₂O

2942 μ l 20XSSC

500 μ l 1M Tris-Cl 7.5

352 5 M NaOH

100 mg DNA (Sigma D-1626)

Dissolve with shaking, heat at 80 C for 10 min. Then neutralize by adding 588 μ l 2M

Tris-Cl (pH 7.4) sonicate at max setting for 1 minute.

Autoclave and aliquot.

50X Denhardt's Solution (100 ml, store at -20 C)

1 g Ficoll (Sigma F4375)

1 g Polyvinylpyrrolidonesigma (Sigma PVP360)

1 g BSA

Filter solution through 0.45 μ filter unit. Aliquot and freeze.

Sodium phosphate, pH 6.8 (100 ml, store at room temp)

46.3 ml 1M Na₂HPO₄ (14.2 g/100ml)

53.7 ml 1M NaH₂PO₄ (13.8 g/100ml)

DNA probe preparation by random labeling

Use Invitrogen kit, #18187-013.

Denature 25 ng DNA in 5-20 μ l water by boiling then place tube on ice. (Use gel isolated fragment from PCR or restriction digest of plasmid)

Add 2 μ l each of dATP, dGTP and dTTP.

Add 15 μ l random primers buffer mixture.

Add 5 μ l (50 μ Ci) of α ³²P dCTP (3000 Ci/mmol)

Add 1 μ l Klenow enzyme, mix gently and centrifuge briefly.

Incubate at room temperature for 1 hour.

Add 5 μ l stop buffer

Purify labeled DNA using Amersham Probe Quant G-50 Micro Column (#27-5335-01) exactly as described by the manufacturer.

Boil probe, cool on ice, then add 5-25 ng to prehyb.

Incubate blot at 42 C overnight in roller bottle.

Blot work-up

Wash twice in 0.1X SSC, 0.1% SDS for 10 min at 42 C in hyb bottle.

Wrap blot in plastic wrap and expose to film or phosphorimager screen. Do not let the blot dry out if you plan to re-probe it.

Strip blot for re-probing

Bring 1 liter 0.1 % SDS just to boiling.

Remove from heat and add blot. Allow solution to cool. Discard solution in radioactive waste. Repeat. Store blots at -20 C.